Equipment - raised toilet seats and toilet aids

1. Raised toilet seat



- Larger front/rear cutouts provide easier access,
- Moulded plastic brackets keep it firmly in place,
- Easily adjusted to fit on most toilet shapes,
- Fitting and removal of the seat is quick,
- Comes in varying heights -2, 4, or 6 inch.

2. Mowbray toilet frame



- Raised toilet seat with frame attached,
- Slots over the toilet,
- Gripped feet to prevent slipping or movement,
- Useful if there is no space to fit wall-fixed rails,
- Useful if there are plasterboard walls,
- Can be a potential trip hazard.

3. Toilet frame



- Can be used with raised toilet seat,
- Slots over the toilet,
- Gripped feet to prevent slipping or movement,
- Useful if there is no space to fit wall-fixed rails,
- Useful if there are plasterboard walls,
- Can be a potential trip hazard.

4. Wall mounted rail



- Can be placed adjacent to toilet,
- Can be placed horizontal or vertical, Provides support to get up from toilet,
- Require a suitable wall for safe fixing.

5. Drop down toilet support rail



- Can be placed to open side of toilet,
- Provides support to get up from toilet, Folds away suitable in compact space, Require suitable wall for safe fixing.